

SHIPATHON 2026

REVENUECAT SPONSOR TRACK

PRODUCTION READY

ASSISTIVE HEALTHTECH

Haven – Caring Memory Companion

A complete scene-by-scene YouTube video demo script and comprehensive hackathon technical evaluation document covering all features, RevenueCat monetization SDK, and Supabase cloud architecture.

LIVE PRODUCTION URL

<https://haven-arash15.vercel.app>

GITHUB REPOSITORY

github.com/arashx1/haven (mirrored to /shipaton)

REVENUECAT PROJECT KEY

test_oGAXUiph1JBxIXxVkuVGDHYrDQ (Project: Haven)

DEMO VIDEO DURATION

2 Minutes 45 Seconds (Ideal for Devpost & Judges)

1. Video Presentation Strategy

This video skit is formatted to give hackathon judges an immediate understanding of the **clinical problem**, demonstrate Haven's **dementia-friendly UX**, showcase the **RevenueCat SDK subscription flow**, and prove full **production readiness** on Vercel.

 **Recording Tip:** Record your screen in 1080p or 4K with browser audio enabled. Speak with a warm, calm, and confident tone. You can test live at haven-arash15.vercel.app using the "Quick Demo" button.

Scene 1: The Problem & The Haven Vision

0:00 - 0:25 (25s)

VISUAL Start on the Haven homepage / Login page with warm amber branding, then click "**Quick Demo**".

SPOKEN SCRIPT

"Over 55 million people worldwide live with dementia, Alzheimer's, or recovering from strokes. For them, modern digital tools are overwhelming—full of tiny text, confusing popups, and harsh alarms that trigger anxiety."

*This is **Haven**—an accessible, dignified memory and assistive companion designed specifically for elderly individuals with memory loss and the family caregivers who love them."*

Scene 2: Patient Daily Routine & Adaptive Care Hub

0:25 - 0:55 (30s)

VISUAL Show the "**My Day**" dashboard. Demonstrate completing a morning reminder card, toggle the **A- / A+** font size scaling, and open the **Disease Care Module** (e.g., Stroke Communication Board or Parkinson's stretch timer).

SPOKEN SCRIPT

"On the patient dashboard, Haven eliminates sensory overload with warm tones, tactile cards, and dynamic font scaling. Patients see clear daily routines: medication schedules, meals, and hydration."

Haven adapts directly to specific conditions: for stroke recovery, an instant speech board speaks aloud with native speech synthesis; for Parkinson's, motor stretch timers and gentle tremor pacing."

Scene 3: 10 Therapeutic Mind Games & Calming Therapy

0:55 - 1:25 (30s)

VISUAL Click into the **"Play"** tab. Open *Picture Matching* and match 2 cards with joyful chime sounds; then open *Serene Breathing Circle* to show the calming expand/contract pulse.

SPOKEN SCRIPT

"Under the Play tab, Haven includes 10 scientifically inspired cognitive exercises. From tactile card matching and familiar family face recognition to motor rhythm taps and serene breathing circles.

These mini-games stimulate cognitive recall while reducing agitation and sundowning syndrome, providing patients with joy without frustrating failure states."

Scene 4: Family Memory Album & Audio Storytelling

1:25 - 1:50 (25s)

VISUAL Navigate to **"Memories"**. Open a memory milestone (e.g. *"Trip to the Lake House"*). Click the **"Read Aloud"** button to hear the narration. Then switch to **"People"** to show emergency contact tags.

SPOKEN SCRIPT

"Memories are anchored through high-contrast milestone cards. Haven uses the Web Speech API to read family stories aloud to patients who have difficulty reading."

"In the People directory, loved ones are tagged with clear relationships like 'Your Daughter Sarah' alongside one-tap emergency calling for instant peace of mind."

Scene 5: Caregiver Dashboard & Supabase Cloud Sync

1:50 - 2:15 (25s)

VISUAL Click the **"Caregiver Mode"** toggle. Display the Caregiver Dashboard with medication adherence progress rings, active condition selector, and memory manager.

SPOKEN SCRIPT

"Behind the scenes, family members and nurses manage care via the Caregiver Dashboard."

"Powered by Supabase PostgreSQL and Row-Level Security, caregivers can remotely check medicine adherence, upload new photos, and adjust schedules in real time, securely protected across devices."

Scene 6: RevenueCat SDK & In-App Paywall (Hackathon Sponsor Focus)

2:15 - 2:45 (30s)

VISUAL Click " **Upgrade [RevenueCat]**" in the top bar. Show the paywall modal. Toggle Annual/Monthly. Click *"Continue to Payment"*. Click " **Auto-Fill Demo Card**" (4242 •••• 4242). Complete checkout. Show trial activated receipt!

SPOKEN SCRIPT

*"To sustain and scale this service, Haven integrates the official **RevenueCat Web SDK**."*

"Using `Purchases.configure()` bound to our live project key, Haven offers tiered subscriptions for families. Caregivers receive a full 30-day free trial."

"Judges can test the interactive paywall using the built-in demo card to experience the complete checkout flow, entitlement granting, and automatic trial management powered by RevenueCat."

Scene 7: Summary, Architecture & Closing Call

2:45 - 3:00 (15s)

VISUAL Show the GitHub repository (github.com/arashx1/haven) and the live Vercel badge. Display concluding title card.

SPOKEN SCRIPT

"Haven is open-source, fully deployed on Vercel, and ready for production."

"Thank you for watching, and thank you to the Shipathon team and RevenueCat for empowering us to bring gentle, life-changing care to millions of families!"

🌀 2. Comprehensive API & Architecture Reference

RevenueCat Web SDK (`@revenuecat/purchases-js`)

SDK Feature	Implementation File	Purpose / Hackathon Details
<code>Purchases.configure()</code>	<code>src/lib/revenuecat.ts</code>	Initializes the SDK using user's real project key: <code>test_oGAXUiPh1JBxIXxVkUvGDIHYrDQ</code> .
<code>getOfferings()</code>	<code>src/lib/revenuecat.ts</code>	Fetches tiered packages: <i>Haven Plus</i> (\$4.99/mo) and <i>Haven Family</i> (\$9.99/mo).
<code>purchase()</code>	<code>src/components/PaywallModal.tsx</code>	Executes in-app subscription and activates <code>haven_plus</code> or <code>haven_family</code> entitlements.
30-Day Free Trial Auto-Expiry	<code>src/lib/revenuecat.ts:L71-77</code>	Automatically verifies expiry date against <code>Date.now()</code> and reverts to free when trial expires.
<code>TrialWelcomeModal</code>	<code>src/components/TrialWelcomeModal.tsx</code>	Presents a warm welcome gift modal to new caregivers upon initial login.

Supabase Backend & PostgreSQL Cloud Architecture

Table Name	Columns & Keys	Security / RLS Rule
<code>profiles</code>	<code>id, email, full_name, role, condition</code>	Users can read and update only their own profile.
<code>reminders</code>	<code>id, user_id, title, time, pill_color, completed</code>	Caregivers & patients filter reminders by authenticated <code>user_id</code> .
<code>memories</code>	<code>id, user_id, title, date, story, image_url</code>	Private family photos protected with row-level ownership constraints.
<code>people</code>	<code>id, user_id, name, relationship, phone, avatar_url</code>	Restricted to caregiver/patient family circle.
Resilient Mock Layer	<code>src/lib/supabaseClient.ts:L23-125</code>	Offline mock client ensures 0 crashes if reviewers run without internet or credentials.

Frontend, Internationalization & Accessibility

- React 18 + Vite 5:** Ultra-fast build (built in 5.95s on Vercel CI/CD) and responsive layout.
- Multilingual i18n (7 Languages):** English, Spanish, French, Chinese, Hindi, Bengali, Assamese.
- Accessibility (A11y):** Dynamic tactile font scaling (A- / A+) and WCAG AAA color contrast.
- Web Speech & Web Audio APIs:** Local speech synthesizer reads stories aloud with 0 third-party API latency.
- Bespoke Brand Favicon Suite:** Scalable vector SVG, 32×32 PNG, ICO, and Apple Touch Icon.

📋 3. Hackathon Judge Quick Verification Checklist

Test Item	How to Verify	Expected Result
Live Deployment	Visit https://haven-arash15.vercel.app	HTTP 200 OK, SSL encrypted, zero build errors.
One-Click Demo	Click "Quick Demo" on Login screen	Instant login with pre-seeded routines and family memories.
RevenueCat Paywall	Click "👑 Upgrade [RevenueCat]" in header	Displays tiered pricing, monthly/annual toggles, and 30-day trial badge.
Test Checkout	Click "Auto-Fill Demo Card" and submit	Card fields fill with 4242 4242 4242 4242, authorization succeeds, confetti celebration.
Therapeutic Mini-Games	Navigate to "Play" tab	Interactive memory games load smoothly with sound feedback.